


Product catalog

Medical Hosiery
Support Hosiery



ORTHOSERVICE
RO+TEN

Take care feel better



Hands typing on a laptop keyboard.



M.S.F.

“

Although we consider our encounter a happy and lucky event, we founded and manage our group willfully. We pursue our objectives everyday with trust, perseverance, and competence.

”

*Pamela Pirovano and Paolo Rossi,
Founders, Swiss Headquarters*

Content

Mission and vision	4
Who we are.....	6
Company system	8

MEDICAL HOSIERY

Natural Line	10
Electa Line	13
Electa Line Kit post op	15
Electa Line Orthopaedic	16
TH Line	19

SUPPORT HOSIERY

Metropolitan Milk Line	22
Suite Lady Line.....	24
Suite Man Line.....	27
Dia Soft Line.....	29

General terms and condition of sales	32
--	----

Every day we begin
to be what we
want to become.





Mission

We design, manufacture and sell, a comprehensive range of orthoses, plantars and products for elastic compression.

Our duo of quality and price, management, and flexibility of production, have positioned us as market leader although we are still a family at heart. Success is not only an individual target but a common achievement of the group.

Vision

We rely on our organization. Our aim is to make our products available to the greatest number of people possible by focusing on flexibility and commitment, in order to be a reference point for orthopedic technicians, experts and patients. Because effective and noninvasive rehabilitation is the means to a better life.

Take
care
feel
better



Italian Design

Designing a product implies getting involved, and everytime it's a new challenge that we accept enthusiastically. We work to make the brace functional as doctors recommend us, practical for technicians that will adapt it and effective in making the patient's life easier. Matching all requirements is not easy, but when we are able to generate design with an original solution which helps a person to live better, it makes it worthwhile for us.

*Francesco Turconi, engineer, R&D,
Swiss Headquarters*



European Manufacturing

Ortoveneta

Here in Italy we are in charge of products having both aluminum and plastic components. However, before applying the components to the products, we have to design machinery for their applications. We need to evoke our creativity by involving our skills such as inventiveness and strictness.

Michele Caporello, CEO, Ortoveneta, Italy

ICS Orthoservice

In the Moldovan craftsmanship we have found a great partner. In Moldavia we manufacture a lot of products that require the utmost accuracy and precision together with quality. Each project is a real chisel work, like true perfectionists.

Maurizio Pilastro, Head of ICS Production, Moldova



Swiss Precision Planning

Putting theory into practice requires intuition and patience.

After checking the project's feasibility, we select raw materials matching our quality standards through preferred European suppliers. We choose the manufacturing site according to the product. Afterwards, we plan terms and timing according to the range of products in our catalogue. Our goal is to define who does what.

*Mauro Ravese, production planning,
Swiss Headquarters*



Market

Orthoservice operates in twenty-nine countries throughout Europe and around the world. Its partners are all the leading distributors to provide guaranteed experience and Swiss quality, prompt support, assistance and product personalization, as we are original equipment manufacturers (OEM). Our headquarters in Switzerland works very closely with the production, research and development, and marketing departments. Opportunities and solutions are promptly assessed with expertise and professionalism, using justified needs and local businesses for the final marketing of the product, ensuring that the brand always maintains its recognizable standard.

*Fabio Albanese, Export and Domestic
Sales Manager*

”



A single company system, many production and management sites



ORTHOERVICE POLSKA SP. Z O.O.
2021 - Częstochowa, Poland



ORTHOERVICE DEUTSCHLAND GmbH
2011 - Baden-Baden, Germany



ICS ORTHOSERVICE
2010 - Balti, Moldova



ORTHOERVICE AG
1996 - Chiasso, Switzerland



ORTOVENETA
2007 - Maser, Italy



RO+TEN S.r.l
1978 - Verano Brianza, Italy

ORTHOERVICE RO+TEN GROUP

- Since 1978 in the orthoses sector
- Over 130 employees
- 9 accredited patents
- Export to 20 countries

Medical hosiery

Graduated compression stockings, socks and tights are surgical medical aids prescribed by a doctor after assessing the needs of an individual patient. By compressing superficial veins and helping blood to return to the deep veins, which usually function normally, they prevent varicose veins from worsening and the onset of problems associated with venous stasis or venous ulcers.

The compression level is highest around the ankle and lessens towards the top of the thigh.

The degree of compression around the ankle is commonly expressed in mmHg (millimetres of mercury).

Orthoservice offers **Class 1** and **Class 2** medical hosiery: they allow blood to return back up from the feet to the heart and from superficial veins to deep veins.



Value of compression at the ankle

- Class 1: 18-21 mmHg
- Class 2: 23-32 mmHg

NATURAL LINE



Unisex

NATURAL



SINGLE PATIENT
MULTIPLE USE



AD Knee Highs
open toe

Cotton blend medical hosiery

Natural graduated compression hosiery is made with the finest quality yarns. The cotton lends the hosiery a natural finish that feels wonderfully soft against the skin. Every product is made with computerised machinery and by hand and has to pass many quality control tests. Comfortable to wear, they retain their compression and do not ladder easily.

INDICATIONS

The compression class is at the discretion of a doctor

- Functional disorders, such as heavy, tired legs, without visible varicose veins
- Mild varicose veins associated with subjective disorders (spider veins, reticular veins)
- Varicose veins during pregnancy with or without oedema
- Localised oedema in the evening (venous)
- Primary varicosities with oedema
- Advanced primary varicosities with oedema (main vein and branch varicosities)
- Chronic venous insufficiency (CVI)
- Prophylaxis of thrombosis and embolisms
- Thrombophlebitis or venous thrombosis (superficial thrombosis) possibly with oedema
- DVT (deep vein thrombosis)
- Post-thrombotic oedema with or without secondary follow-up therapy
- To ensure successful therapy:
 1. following sclerotherapy
 2. following an elective outpatient procedure for various veins
 3. following stripping
- Before sclerotherapy or surgery to prevent the onset of oedema
- Ulcer prevention and therapy
- Postoperativ oedema e
- Post-traumatic oedema after basic therapy
- Reversible lymphoedema following complex decongestive therapy
- Irreversible lymphoedema with skin becoming hard and tight following complex decongestive therapy
- Lipoedema of grade II and above
- Idiopathic cyclic oedema
- Lymphoedema following basic therapy
- Congestive disorders caused by immobility
- Angiodysplasia

CONTRAINDICATIONS

- Advanced peripheral arterial occlusive disease
- Decompensated heart failure
- Septic phlebitis (venous inflammation caused by microorganisms)
- Phlegmasia cerulea dolens
- Weeping dermatosis
- Intolerance to the material used for graduated compression hosiery
- Loss of sensation in the extremities
- Severe peripheral neuropathy (e.g. diabetes mellitus)
- Chronic polyarthritis

The company accepts no responsibility for improper use of a product

FEATURES

- Cotton blend
- Medical hosiery made on new-generation computerised circular knitting machines
- Any seams are stretchy and extremely fine
- Made with the finest quality Swiss yarns:
 - fine
 - soft
 - guaranteed to be long-lasting and resistant to ointments and creams
 - excellent for maintaining compression and elasticity over time
- Yarn-dyed for more even colours and better long-term elasticity
- A special weave is used to make the toe and heel even stronger
- Extensive quality control involves visual inspection and testing by equipment with sensors



AG One-leg Tights
open toe



AT Tights



AG Stockings
open toe

COMPOSITION

- Class 1: 50% Polyamide - 30% Cotton - 20% Elastane
- Class 2: 46% Polyamide - 30% Cotton - 24% Elastane

Value of compression at the ankle

- Class 1: 18-21 mmHg
- Class 2: 23-32 mmHg

AVAILABLE STYLES

NATURAL LINE	Compression class		Colour	Toe		Length	
	1	2	Beige	Open	Closed	Standard	Long
AD Knee Highs	•	•	•	•	ND	•	•
AG Stockings	•	•	•	•	ND	•	•
AG One-leg Tights (right/left)	•	•	•	•	ND	•	•
AT Tights	•	•	•	ND	•	•	•


SIZE GUIDE

CIRCUMFERENCE (cm)

SIZE		XS	S	M	L	XL	XXL
Waist	cT	60-70	70-80	80-90	90-100	100-110	110-120
Thigh	cG	48-56	52-60	56-64	60-68	64-72	68-76
Mid Thigh	cF	39-44	43-48	47-52	51-56	55-60	59-64
Below the knee	cD	26-31	30-35	34-39	38-43	42-47	46-51
Calf	cC	29-34	32-37	35-40	38-43	42-47	46-51
Ankle	cB	19-21	21-23	23-25	25-27	27-29	29-31

LENGTH (CM)

STYLE	SIZE	STANDARD	LONG
AD Knee Highs	/D	35-41	40-45
AG Stockings	/G	65-72	73-82
AG One-leg Tights	/G	65-72	73-82
AT Tights	/G	65-72	73-82



LINE	Code	Size	Range	Toe	Lenght	Compression	Colour
Natural	MNA059	XS-S-M-L-XL-XXL	AD Knee Highs	open	Standard	1 CCL	Beige
Natural	MNA060	XS-S-M-L-XL-XXL	AD Knee Highs	open	Long	1 CCL	Beige
Natural	MNA061	XS-S-M-L-XL-XXL	AG Stockings	open	Standard	1 CCL	Beige
Natural	MNA062	XS-S-M-L-XL-XXL	AG Stockings	open	Long	1 CCL	Beige
Natural	MNA098	XS-S-M-L-XL-XXL	AG One-leg Tights (right)	open	Standard	1 CCL	Beige
Natural	MNA100	XS-S-M-L-XL-XXL	AG One-leg Tights (left)	open	Standard	1 CCL	Beige
Natural	MNA099	XS-S-M-L-XL-XXL	AG One-leg Tights (right)	open	Long	1 CCL	Beige
Natural	MNA101	XS-S-M-L-XL-XXL	AG One-leg Tights (left)	open	Long	1 CCL	Beige
Natural	MNA106	XS-S-M-L-XL-XXL	AT Tights	closed	Standard	1 CCL	Beige
Natural	MNA107	XS-S-M-L-XL-XXL	AT Tights	closed	Long	1 CCL	Beige
Natural	MNA063	XS-S-M-L-XL-XXL	AD Knee Highs	open	Standard	2 CCL	Beige
Natural	MNA064	XS-S-M-L-XL-XXL	AD Knee Highs	open	Long	2 CCL	Beige
Natural	MNA065	XS-S-M-L-XL-XXL	AG Stockings	open	Standard	2 CCL	Beige
Natural	MNA066	XS-S-M-L-XL-XXL	AG Stockings	open	Long	2 CCL	Beige
Natural	MNA102	XS-S-M-L-XL-XXL	AG One-leg Tights (right)	open	Standard	2 CCL	Beige
Natural	MNA104	XS-S-M-L-XL-XXL	AG One-leg Tights (left)	open	Standard	2 CCL	Beige
Natural	MNA103	XS-S-M-L-XL-XXL	AG One-leg Tights (right)	open	Long	2 CCL	Beige
Natural	MNA105	XS-S-M-L-XL-XXL	AG One-leg Tights (left)	open	Long	2 CCL	Beige
Natural	MNA106	XS-S-M-L-XL-XXL	AT Tights	closed	Standard	2 CCL	Beige
Natural	MNA109	XS-S-M-L-XL-XXL	AT Tights	closed	Long	2 CCL	Beige

ELECTA LINE



Unisex

ELECTA



AT Tights

FEATURES

- Made with multi-fibre yarn
- Medical hosiery made on new-generation computerised circular knitting machines
- Any seams are stretchy and extremely fine
- Made with the finest quality Swiss yarns:
 - fine
 - soft
 - guaranteed to be long-lasting and resistant to ointments and creams
- - excellent for maintaining compression and elasticity over time
- Yarn-dyed for more even colours and better long-term elasticity
- A special weave is used to make the toe and heel even stronger
- Extensive quality control involves visual inspection and testing by equipment with sensors

Multi-fibre yarn medical hosiery



SINGLE PATIENT
MULTIPLE USE

Electa graduated compression hosiery blends successful therapy with elegance and style. Soft and fine, they ensure women can always dress with class on any occasion.

Every product is made with computerised machinery and by hand and has to pass many quality control tests.

Comfortable to wear, they retain their compression and do not ladder easily.

INDICATIONS

The compression class is at the discretion of a doctor

- Functional disorders, such as heavy, tired legs, without visible varicose veins
- Mild varicose veins associated with subjective disorders (spider veins, reticular veins)
- Varicose veins during pregnancy with or without oedema
- Localised oedema in the evening (venous)
- Primary varicosities with oedema
- Advanced primary varicosities with oedema (main vein and branch varicosities)
- Chronic venous insufficiency (CVI)
- Prophylaxis of thrombosis and embolisms
- Thrombophlebitis or venous thrombosis (superficial thrombosis) possibly with oedema
- DVT (deep vein thrombosis)
- Post-thrombotic oedema with or without secondary follow-up therapy
- To ensure successful therapy:
 1. following sclerotherapy
 2. following an elective outpatient procedure for various veins
 3. following stripping
- Before sclerotherapy or surgery to prevent the onset of oedema
- Ulcer prevention and therapy
- Postoperative oedema
- Post-traumatic oedema after basic therapy
- Reversible lymphoedema following complex decongestive therapy
- Irreversible lymphoedema with skin becoming hard and tight following complex decongestive therapy
- Lipoedema of grade II and above
- Idiopathic cyclic oedema
- Lymphoedema following basic therapy
- Congestive disorders caused by immobility
- Angiodysplasi


CONTRAINDICATIONS

- Advanced peripheral arterial occlusive disease
- Decompensated heart failure
- Septic phlebitis (venous inflammation caused by microorganisms)
- Phlegmasia cerulea dolens
- Weeping dermatosis
- Intolerance to the material used for graduated compression hosiery
- Loss of sensation in the extremities
- Severe peripheral neuropathy (e.g. diabetes mellitus)
- Chronic polyarthritis

The company accepts no responsibility for improper use of a product.



AG Stockings
open toe



AG One-leg Tights
open toe



AD Knee Highs
closed toe

COMPOSITION

- Class 1: 74% Polyamide - 26% Elastane
- Class 2: 70% Polyamide - 30% Elastane

Value of compression at the ankle:

- Class 1: 18-21 mmHg
- Class 2: 23-32 mmHg

AVAILABLE STYLES

ELECTA LINE	Colour	Compression class		Toe		Length	
		1	2	Open	Closed	Standard	Long
AD Knee Highs	Beige Black	• •	• •	• •	• •	• •	• •
AG Stockings	Beige Black	• •	• •	• •	• (only CL1) •	• •	• •
AG One-leg Tights (right/left)	Beige	•	•	•	ND	•	•
AT Tights	Beige Black	• •	• •	ND ND	• •	• •	• •


SIZE GUIDE

CIRCUMFERENCE (cm)

SIZE		XS	S	M	L	XL	XXL
Waist	cT	60-70	70-80	80-90	90-100	100-110	110-120
Thigh	cG	48-56	52-60	56-64	60-68	64-72	68-76
Mid Thigh	cF	39-44	43-48	47-52	51-56	55-60	59-64
Below the knee	cD	26-31	30-35	34-39	38-43	42-47	46-51
Calf	cC	29-34	32-37	35-40	38-43	42-47	46-51
Ankle	cB	19-21	21-23	23-25	25-27	27-29	29-31

LENGTH (CM)

STYLE	SIZE	STANDARD	LONG
AD Knee Highs	/D	35-41	40-45
AG Stockings	/G	65-72	73-82
AG One-leg Tights	/G	65-72	73-82
AT Tights	/G	65-72	73-82



ELECTA LINE



Unisex

ELECTA



SINGLE PATIENT
MULTIPLE USE

Kit post op

This kit contains Electa one-leg tights and two liners and is specially designed for the delicate postoperative period following surgery on a lower limb. They provide greater compression around than thigh than traditional one-leg tights which helps to reduce any oedema in the upper part of the leg. The outstanding elasticity of the multi-fibre yarn, special weave and excellent breathability combine guaranteed efficacy and maximum possible comfort to make this the ideal hosiery to wear during the day and at night. The fine liners are gentle on the skin and help keep dressings firmly in place, making it even easier to put on the tights correctly.

AVAILABLE STYLES

ELECTA LINE POST OP	Colour	Compression class		Toe		Lenght	
		1	2	Open	Closed	Standard	Long
AG One-leg Tights (right/left)	Beige	•	•	•	ND	•	•

SIZE GUIDE

CIRCUMFERENCE (cm)

SIZE		XS	S	M	L	XL	XXL
Waist	cT	60-70	70-80	80-90	90-100	100-110	110-120
Thigh	cG	48-56	52-60	56-64	60-68	64-72	68-76
Mid Thigh	cF	39-44	43-48	47-52	51-56	55-60	59-64
Below the knee	cD	26-31	30-35	34-39	38-43	42-47	46-51
Calf	cC	29-34	32-37	35-40	38-43	42-47	46-51
Ankle	cB	19-21	21-23	23-25	25-27	27-29	29-31

LENGTH (CM)

STYLE	SIZE	STANDARD	LONG
AG One-leg Tights	/G	65-72	73-82

COMPOSITION

- Class 1: 50% Polyamide 30% Cotton - 20% Elastane
- Class 2: 46% Polyamide 30% Cotton - 24% Elastane

Value of compression at the ankle:

- Class 1: 18-21 mmHg
- Class 2: 23-32 mmHg

Skin-friendly understocking to keep in place the dressing

AG One-leg Tights

CONTENTS:

AG One-leg Tights CL2 open toe (standard o long)
2 liners for dressings (one size)




ELECTA LINE



Unisex

ELECTA

AD Knee Highs
open toe



Electa Orthopaedic

The single knee-high is ideal for patients with oedema following surgery, the removal of plaster of Paris or a trauma.



SINGLE PATIENT
MULTIPLE USE

AVAILABLE STYLES

ELECTA LINE ORTHOPAEDIC	Compression class		Colour	Toe		Length	
	1	2	Beige	Open	Closed	Standard	Long
AD Knee Highs	•	•	•	•	ND	•	•

SIZE GUIDE

CIRCUMFERENCE (cm)

SIZE		XS	S	M	L	XL	XXL
Below the knee	cD	26-31	30-35	34-39	38-43	42-47	46-51
Calf	cC	29-34	32-37	35-40	38-43	42-47	46-51
Ankle	cB	19-21	21-23	23-25	25-27	27-29	29-31

LENGTH (CM)


STYLE	SIZE	STANDARD	LONG
AD Knee Highs	/D	35-41	40-45

COMPOSITION

- Class 1: 74% Polyamide - 26% Elastane
- Class 2: 70% Polyamide - 30% Elastane

Value of compression at the ankle:

- Class 1: 18-21 mmHg
- Class 2: 23-32 mmHg



LINE	Code	Size	Range	Toe	Lenght	Compression	Colour
Electa	MEL067	XS-S-M-L-XL-XXL	AD Knee Highs	open	Standard	1 CCL	Beige
Electa	MEL128	XS-S-M-L-XL-XXL	AD Knee Highs	open	Standard	1 CCL	Black
Electa	MEL068	XS-S-M-L-XL-XXL	AD Knee Highs	open	Long	1 CCL	Beige
Electa	MEL129	XS-S-M-L-XL-XXL	AD Knee Highs	open	Long	1 CCL	Black
Electa	MEL069	XS-S-M-L-XL-XXL	AD Knee Highs	closed	Standard	1 CCL	Beige
Electa	MEL112	XS-S-M-L-XL-XXL	AD Knee Highs	closed	Standard	1 CCL	Black
Electa	MEL070	XS-S-M-L-XL-XXL	AD Knee Highs	closed	Long	1 CCL	Beige
Electa	MEL113	XS-S-M-L-XL-XXL	AD Knee Highs	closed	Long	1 CCL	Black
Electa	MEL075	XS-S-M-L-XL-XXL	AG Stockings	open	Standard	1 CCL	Beige
Electa	MEL132	XS-S-M-L-XL-XXL	AG Stockings	open	Standard	1 CCL	Black
Electa	MEL076	XS-S-M-L-XL-XXL	AG Stockings	open	Long	1 CCL	Beige
Electa	MEL133	XS-S-M-L-XL-XXL	AG Stockings	open	Long	1 CCL	Black
Electa	MEL077	XS-S-M-L-XL-XXL	AG Stockings	closed	Standard	1 CCL	Beige
Electa	MEL116	XS-S-M-L-XL-XXL	AG Stockings	closed	Standard	1 CCL	Black
Electa	MEL078	XS-S-M-L-XL-XXL	AG Stockings	closed	Long	1 CCL	Beige
Electa	MEL117	XS-S-M-L-XL-XXL	AG Stockings	closed	Long	1 CCL	Black
Electa	MEL081	XS-S-M-L-XL-XXL	AG One-leg Tights (right)	open	Standard	1 CCL	Beige
Electa	MEL082	XS-S-M-L-XL-XXL	AG One-leg Tights (right)	open	Long	1 CCL	Beige
Electa	MEL083	XS-S-M-L-XL-XXL	AG One-leg Tights (left)	open	Standard	1 CCL	Beige
Electa	MEL084	XS-S-M-L-XL-XXL	AG One-leg Tights (left)	open	Long	1 CCL	Beige
Electa	MEL089	XS-S-M-L-XL-XXL	AT Tights	closed	Standard	1 CCL	Beige
Electa	MEL120	XS-S-M-L-XL-XXL	AT Tights	closed	Standard	1 CCL	Black
Electa	MEL090	XS-S-M-L-XL-XXL	AT Tights	closed	Long	1 CCL	Beige
Electa	MEL121	XS-S-M-L-XL-XXL	AT Tights	closed	Long	1 CCL	Black
Electa	MEL071	XS-S-M-L-XL-XXL	AD Knee Highs	open	Standard	2 CCL	Beige
Electa	MEL130	XS-S-M-L-XL-XXL	AD Knee Highs	open	Standard	2 CCL	Black
Electa	MEL072	XS-S-M-L-XL-XXL	AD Knee Highs	open	Long	2 CCL	Beige
Electa	MEL131	XS-S-M-L-XL-XXL	AD Knee Highs	open	Long	2 CCL	Black
Electa	MEL073	XS-S-M-L-XL-XXL	AD Knee Highs	closed	Standard	2 CCL	Beige
Electa	MEL114	XS-S-M-L-XL-XXL	AD Knee Highs	closed	Standard	2 CCL	Black
Electa	MEL074	XS-S-M-L-XL-XXL	AD Knee Highs	closed	Long	2 CCL	Beige
Electa	MEL115	XS-S-M-L-XL-XXL	AD Knee Highs	closed	Long	2 CCL	Black
Electa	MEL079	XS-S-M-L-XL-XXL	AG Stockings	open	Standard	2 CCL	Beige
Electa	MEL134	XS-S-M-L-XL-XXL	AG Stockings	open	Standard	2 CCL	Black
Electa	MEL080	XS-S-M-L-XL-XXL	AG Stockings	open	Long	2 CCL	Beige
Electa	MEL135	XS-S-M-L-XL-XXL	AG Stockings	open	Long	2 CCL	Black
Electa	MEL118	XS-S-M-L-XL-XXL	AG Stockings	closed	Standard	2 CCL	Black
Electa	MEL119	XS-S-M-L-XL-XXL	AG Stockings	closed	Long	2 CCL	Black
Electa	MEL085	XS-S-M-L-XL-XXL	AG One-leg Tights (right)	open	Standard	2 CCL	Beige
Electa	MEL086	XS-S-M-L-XL-XXL	AG One-leg Tights (right)	open	Long	2 CCL	Beige
Electa	MEL087	XS-S-M-L-XL-XXL	AG One-leg Tights (left)	open	Standard	2 CCL	Beige
Electa	MEL088	XS-S-M-L-XL-XXL	AG One-leg Tights (left)	open	Long	2 CCL	Beige
Electa	MEL093	XS-S-M-L-XL-XXL	AT Tights	closed	Standard	2 CCL	Beige
Electa	MEL122	XS-S-M-L-XL-XXL	AT Tights	closed	Standard	2 CCL	Black
Electa	MEL094	XS-S-M-L-XL-XXL	AT Tights	closed	Long	2 CCL	Beige
Electa	MEL123	XS-S-M-L-XL-XXL	AT Tights	closed	Long	2 CCL	Black

LINE	Code	Size	Range	Toe	Lenght	Compression	Colour
Electa Orthopaedic	MER124	XS-S-M-L-XL-XXL	AD Knee Highs	open	Standard	1 CCL	Beige
Electa Orthopaedic	MER125	XS-S-M-L-XL-XXL	AD Knee Highs	open	Long	1 CCL	Beige
Electa Orthopaedic	MER126	XS-S-M-L-XL-XXL	AD Knee Highs	open	Standard	2 CCL	Beige
Electa Orthopaedic	MER127	XS-S-M-L-XL-XXL	AD Knee Highs	open	Long	2 CCL	Beige
Electa Post Op	MEK136	XS-S-M-L-XL-XXL	AG One-leg Tights (right)	open	Standard	2 CCL	Beige
Electa Post Op	MEK138	XS-S-M-L-XL-XXL	AG One-leg Tights (left)	open	Standard	2 CCL	Beige
Electa Post Op	MEK137	XS-S-M-L-XL-XXL	AG One-leg Tights (right)	open	Long	2 CCL	Beige
Electa Post Op	MEK139	XS-S-M-L-XL-XXL	AG One-leg Tights (left)	open	Long	2 CCL	Beige

TH LINE



Unisex



SINGLE PATIENT
MULTIPLE USE



LATEX FREE

Anti-embolism stockings

MEDICAL HOSIERY • TH LINE



AG One-leg Tights Openable



AD Knee Highs open toe



AG One-leg Tights

COMPOSITION

- 82% Polyamide - 18% Elastane

AVAILABLE STYLES

TH LINE	Colour	Toe		Length	
	White	Open	Openable	Standard	Long
AD Knee Highs	•	•	ND	•	•
AG Stockings	•	•	ND	•	•
AG One-leg Tights (right/left)	•	ND	•	•	•

SIZE GUIDE

KNEE HIGH SIZE GUIDE

SIZE	S	M	L	XL
cB [cm]	18-21	21-24	24-27	27-30
cC [cm]	26-32	32-38	38-44	44-50
ID [cm]	Standard: up to 41			
	Long: over 41			


STOCKINGS SIZE GUIDE

SIZE	S	M	L	XL
cB [cm]	18-21	21-24	24-27	27-30
cC [cm]	26-32	32-38	38-44	44-50
cF [cm]	45-55	50-60	55-65	60-70
IG [cm]	Standard: up to 74			
	Long: over 74			

ONE-LEG TIGHTS SIZE GUIDE*

SIZE	S	M	L	XL
cB [cm]	18-21	21-24	24-27	27-30
cC [cm]	26-32	32-38	38-44	44-50
cF [cm]	45-55	50-60	55-65	60-70
IG [cm]	Standard: up to 74			
	Long: over 74			

* Specify right or left



LINE	Code	Size	Range	Toe	Lenght	Compression	Colour
TH	MTH051	S-M-L-XL	AD Knee Highs	open	Standard	18/23 mmHg	White
TH	MTH052	S-M-L-XL	AD Knee Highs	open	Long	18/23 mmHg	White
TH	MTH053	S-M-L-XL	AG Stockings	open	Standard	18/23 mmHg	White
TH	MTH054	S-M-L-XL	AG Stockings	open	Long	18/23 mmHg	White
TH	MTH055	S-M-L-XL	AG One-leg Tights (right)	openable	Standard	18/23 mmHg	White
TH	MTH056	S-M-L-XL	AG One-leg Tights (right)	openable	Long	18/23 mmHg	White
TH	MTH057	S-M-L-XL	AG One-leg Tights (left)	openable	Standard	18/23 mmHg	White
TH	MTH058	S-M-L-XL	AG One-leg Tights (left)	openable	Long	18/23 mmHg	White

Support hosiery

Compression is the most effective conservative **therapy** for offering relief when legs feel tired and heavy, which often occurs after standing for long periods of time. The muscles in the lower legs “naturally pump” blood back up towards the heart. When muscles are inactive for extended periods blood has to beat gravity with circulation alone; at this point venous pressure can increase up to ten times normal levels and in the long term the valves that prevent downward reflux become fatigued. Wearing compression hosiery is an easy solution to implement and is extremely helpful for keeping legs feeling toned, supple and healthy.

METROPOLITAN MILK LINE



Unisex



Milk fibre yarn support hosiery



SINGLE PATIENT
MULTIPLE USE

The Metropolitan Milk collection of knee-high socks is specially designed for people with particularly sensitive skin.

Made from a special yarn featuring a blend of viscose and milk protein fibre (natural fibre), they feel exceptionally soft against the skin yet still provide all the compression you need.

Thanks to the natural fibres, the hosiery absorbs sweat, is extremely comfortable and completely compatible with all skin types.

Anatomically designed reinforced heel and toe. The comfort-fit tops do not restrict circulation.

INDICATIONS

Ideal for

- mild chronic venous insufficiency,
- varicose veins,
- swelling

FEATURES

- The natural milk protein fibre feels soft and is compatible with all skin types.
- Anatomically designed reinforced heel and toe.
- Extensive quality control involves visual inspection and testing by equipment with sensors

COMPOSITION

- 53% Viscose - 23% (Milk) protein - 14% Polyamide - 10% Elastane

AVAILABLE STYLES							
METROPOLITAN MILK LINE	Compression	Colour		Toe		Length	
	15/18 mmHg	Black	Blue	Open	Closed	Standard	Long
Knee Highs	•	•	•	ND	•	•	ND

SIZE GUIDE					
SIZE	2/S	3/M	4/L	5/XL	6/XXL
Shoe size [EU]	37-39	39-41	41-43	43-45	45-47
Foot measurement [cm]	23-25	24.4-26.4	25.7-27.7	27-29	28.4-30.4
Ankle circumference [cm]	20-22	22-24	24-26	26-28	28-30

Value of compression at the ankle:

- 15-18 mmHg

LINE	Code	Size	Range	Toe	Lenght	Compression	Colour
Metropolitan Milk	MMM015	S-M-L-XL-XXL	Knee Highs	closed	Standard	15/18 mmHg	Black
Metropolitan Milk	MMM016	S-M-L-XL-XXL	Knee Highs	closed	Standard	15/18 mmHg	Blue



SUITE LADY LINE



Lady

SUITE
LADY



Tights

Mesh knit support hosier



SINGLE PATIENT
MULTIPLE USE

The Suite Lady hosiery collection is designed for ladies seeking a discreet way to keep their legs healthy.

Just as sheer, elegant and comfortable as any fashionable tights or stockings, they offer the healthy added benefit of compression.

The veins and lymphatic system are supported throughout the day whether you are walking, standing or sitting, giving you an immediate feeling of relief and preventing swelling in the evening. Anatomically designed reinforced heel and toe.

INDICATIONS

12-14mmHg for:

- minor circulation problems
- a tendency to suffer from varicose veins
- during pregnancy

15-18mmHg for:

- mild chronic venous insufficiency
- varicose veins
- swelling
- during pregnancy

FEATURES

- Mesh knit.
- Anatomically designed reinforced heel and toe.
- Extensive quality control involves visual inspection and testing by equipment with sensors

COMPOSITION

- 12-14mmHg (70 DEN): 83% Polyamide, 17% Elastane
- 15-18mmHg (140 DEN): 75% Polyamide, 25% Elastane



Stockings



Knee highs

AVAILABLE STYLES

SUITE LADY LINE	Compression class		Colour			Toe		Length	
	12/14mmHg (70 DEN)	15/18mmHg (140 DEN)	Black	Playa	Clear	Open	Closed	Standard	Long
Knee highs	•	•	•	•	•	ND	•	•	ND
Stockings	•	•	•	•	•	ND	•	•	ND
Tights	•	•	•	•	•	ND	•	•	ND
Tights with boxer brief	•	ND	ND	ND	•	ND	•	•	ND
	ND	•	•	ND	ND	ND	•	•	ND

SIZE GUIDE

CIRCUMFERENCE (cm)

SIZE	1	2	3	4	5
Height [cm]	154-162	158-166	162-174	166-182	174-190
Weight [kg]	43-53	48-63	58-73	63-83	68-93
Ankle circumference [cm]	18-20.5	20.5-23	23-25.5	25.5-28	28-30.5



LINE	Code	Size	Range	Toe	Compression	Colour
Suite Lady	MSL017	1-2-3-4-5	Knee Highs	closed	12/14 mmHg (70 DEN)	Black
Suite Lady	MSL026	1-2-3-4-5	Knee Highs	closed	12/14 mmHg (70 DEN)	Playa
Suite Lady	MSL023	1-2-3-4-5	Knee Highs	closed	12/14 mmHg (70 DEN)	Clear
Suite Lady	MSL018	1-2-3-4-5	Stockings	closed	12/14 mmHg (70 DEN)	Black
Suite Lady	MSL027	1-2-3-4-5	Stockings	closed	12/14 mmHg (70 DEN)	Playa
Suite Lady	MSL024	1-2-3-4-5	Stockings	closed	12/14 mmHg (70 DEN)	Clear
Suite Lady	MSL019	1-2-3-4-5	Tights	closed	12/14 mmHg (70 DEN)	Black
Suite Lady	MSL028	1-2-3-4-5	Tights	closed	12/14 mmHg (70 DEN)	Playa
Suite Lady	MSL025	1-2-3-4-5	Tights	closed	12/14 mmHg (70 DEN)	Clear
Suite Lady	MSL031	1-2-3-4-5	Knee Highs	closed	15/18 mmHg (140 DEN)	Black
Suite Lady	MSL040	1-2-3-4-5	Knee Highs	closed	15/18 mmHg (140 DEN)	Playa
Suite Lady	MSL037	1-2-3-4-5	Knee Highs	closed	15/18 mmHg (140 DEN)	Clear
Suite Lady	MSL032	1-2-3-4-5	Stockings	closed	15/18 mmHg (140 DEN)	Black
Suite Lady	MSL041	1-2-3-4-5	Stockings	closed	15/18 mmHg (140 DEN)	Playa
Suite Lady	MSL038	1-2-3-4-5	Stockings	closed	15/18 mmHg (140 DEN)	Clear
Suite Lady	MSL033	1-2-3-4-5	Tights	closed	15/18 mmHg (140 DEN)	Black
Suite Lady	MSL042	1-2-3-4-5	Tights	closed	15/18 mmHg (140 DEN)	Playa
Suite Lady	MSL039	1-2-3-4-5	Tights	closed	15/18 mmHg (140 DEN)	Clear

SUITE MAN LINE




Man



SINGLE PATIENT
MULTIPLE USE

SUPPORT HOSIERY • SUITE MAN LINE

SUITE
MAN



Cotton blend support socks

The Suite Man collection of knee-high socks is designed for today's modern men who are often sitting or on their feet for hours at a time, or on long-haul flights and coping with the discomfort of pressurisation and sitting in a confined space for long periods.

Made with 70% cotton, the knee-high socks feel incredibly soft against the skin yet still provided all the beneficial compression you need.

Anatomically designed reinforced heel and toe. The comfort-fit tops do not restrict circulation.

INDICATIONS

12-14mmHg for:

- minor circulation problems
- a tendency to suffer from varicose veins;

15-18mmHg for:

- mild chronic venous insufficiency
- varicose veins
- swelling

FEATURES

- 70% cotton.
- Reinforced calf area.
- Anatomically designed reinforced heel and toe.
- Extensive quality control involves visual inspection and testing by equipment with sensors.

COMPOSITION

- 70% Cotton - 15% Polyamide - 15% Elastane



AVAILABLE STYLES

SUITE MAN LINE	Compression class		Colour		Toe		Length	
	12/14mmHg	15/18mmHg	Black	Blue	Open	Closed	Standard	Long
Knee Highs	•	•	•	•	ND	•	•	ND

SIZE GUIDE

SIZE	2/S	3/M	4/L	5/XL	6/XXL
Shoe size [EU]	37-39	39-41	41-43	43-45	45-47
Foot measurement [cm]	23-25	24.4-26.4	25.7-27.7	27-29	28.4-30.4
Ankle circumference [cm]	20-22	22-24	24-26	26-28	28-30

LINE	Code	Size	Range	Toe	Compression	Colour
Suite Man	MSM047	S-M-L-XL-XXL	Knee Highs	closed	12/14 mmHg	Black
Suite Man	MSM048	S-M-L-XL-XXL	Knee Highs	closed	12/14 mmHg	Blue
Suite Man	MSM049	S-M-L-XL-XXL	Knee Highs	closed	15/18 mmHg	Black
Suite Man	MSM050	S-M-L-XL-XXL	Knee Highs	closed	15/18 mmHg	Blue



DIA SOFT LINE



Unisex



SINGLE PATIENT
MULTIPLE USE

SUPPORT HOSIERY • SUITE MAN LINE



AD Knee Highs

Crabyon diabetic socks

INDICATIONS

Ideal for:

- diabetic foot
- arthritis
- athlete's foot
- sensitive feet

FEATURES

- Ultrafine
- Seamless
- Breathable and antibacterial Crabyon functional fibre
- Extensive quality control involves visual inspection and testing by equipment with sensors

COMPOSITION

- 66% Cotton - 17% viscose/Crabyon - 11% Polyamide - 6% Elastane




Mid-calf Socks

AVAILABLE STYLES					
DIA SOFT LINE	Colour			Toe	
	White	Black	Blue	Open	Closed
Mid-calf Socks	•	•	•	ND	•
AD Knee Highs	ND	•	•	ND	•

SIZE GUIDE						
SIZE	1/XS	2/S	3/M	4/L	5/XL	6/XXL
Shoe size [EU]	35-37	37-39	39-41	41-43	43-45	45-47
Foot measurement [cm]	21.7-23.7	23-25	24.4-26.4	25.7-27.7	27-29	28.4-30.4

LINE	Code	Size	Range	Toe	Compression	Colour
DiaSoft	MDS001	XS-S-M-L-XL-XXL	Mid-calf Socks	closed		Black
DiaSoft	MDS002	XS-S-M-L-XL-XXL	Mid-calf Socks	closed		Blue
DiaSoft	MDS095	XS-S-M-L-XL-XXL	Mid-calf Socks	closed		White
DiaSoft	MDS096	XS-S-M-L-XL-XXL	AD Knee Highs	closed		Black
DiaSoft	MDS097	XS-S-M-L-XL-XXL	AD Knee Highs	closed		Blue



Antibacterial, hypoallergenic

The Crabylon textile fiber used for Dia Soft knee-highs is derived from crab chitin. It is antibacterial and antimicrobial, haemostatic, completely biodegradable, hypoallergenic, ecological and biocompatible, these properties remain unchanged even after numerous washes. Crabylon has a high moisture absorbing power and in contact with the skin it is able to prevent dehydration, keeping it fresh and free from bad odors for a relatively long time.

GENERAL TERMS OF SALES

1. Entire Agreement

This document contains all of the terms and conditions relating to the purchase of goods and/or services (collectively, "Goods") by ORTHOSERVICE AG and its affiliates ("Buyer"), from the supplier identified herein ("Seller") and constitutes a legally binding agreement between ORTHOSERVICE AG and Seller ("Agreement"). Buyer and Seller are individually referred to herein as a "Party" and collectively as the "Parties".

2. Orders

From time to time Buyer may submit an order ("Order") for Goods with Seller. Any Order submitted by Buyer to Seller shall be subject to and deemed to incorporate the terms and conditions of this Agreement. Unless agreed to in writing by a duly authorized representative of Buyer, Buyer objects to, and is not bound by, any term or condition that differs from or adds to Buyer's Order or the terms and conditions herein. Commencement of performance by Seller or shipment of any part of the Goods covered by an Order shall be deemed to constitute assent to the Order and the terms and conditions herein.

3. Price

Buyer shall pay Seller the price for the Goods shown in the Order. The price is firm and not subject to increase without Seller's prior written communication. Price will be displayed in the Order Confirmation of Seller for final approval of Buyer.

4. Taxes

Unless specified in an Order, the prices stated do not include sales, use or similar taxes directly applicable to the Goods sold to Buyer. Notwithstanding the foregoing, Buyer shall only be responsible for such taxes when Seller is required by law to collect such taxes from Buyer. Such taxes shall be separately stated on Seller's invoice. Seller shall not add such taxes to Buyer's invoice if Buyer has furnished a valid tax exemption certificate to Seller.

5. Payment

Unless otherwise stated in the Supplier Information document attached hereto or different agreement, Buyer shall prepay Seller the price for the Goods. Seller shall send Buyer a detailed invoice for each Order. Buyer shall make payment to Seller by electronic funds transfer.

6. Delivery

Seller shall deliver Goods in accordance with the schedule in Buyer's Order. Time is of the essence in the performance of the Order. Buyer may select the mode of transportation and the carrier for the Goods. Seller shall pre-pay shipping charges and add such charges as a separate line item on the invoice, unless otherwise agreed in writing. If Buyer requests that Goods be shipped "collect", Seller shall use a Buyer approved transportation carrier and Buyer shall not pre-pay shipping charges. Seller shall be responsible for excess transportation costs resulting from any deviation from Buyer's shipping instructions. If, prior to the delivery of the Goods, Seller has reason to believe that it will be unable to meet the schedule in Buyer's Order, Seller shall immediately notify Buyer in writing, indicate the cause of delay, and use its best efforts to cure the anticipated delay. Upon receipt of notice of the anticipated delay or upon occurrence of an actual delay in performance, Buyer may (i) direct expedited shipment of the Goods with excess costs paid by Seller, or (ii) cancel the Order. Buyer may return at Seller's expense Goods delivered early to Buyer. Goods shall be delivered by Seller to the delivery location specified in the Order or, if delivery location is not specified in the Order, then to Buyer's place of business from which the Goods were ordered. Risk of loss of the Goods while in transit is regulated by the usual international European regulations.

7. Packing

The Goods to be furnished on the Order shall be suitably packed to protect the Goods from damage, prepared for shipment to secure the lowest transportation rates possible (unless a premium shipment method is specified on the Order) and comply with carrier

regulations. No charges shall be allowed for packing, crating, expedited shipping or cartage unless authorized on the Order.

8. Returns

Seller applies merchandise return policies for defective goods only after Seller investigation and written confirmation, which do not comply with European Community standards regulation, and/or for deliveries of products different from the order. In this case only, the cost of transport is belonging to the Seller.

For any other reason for returning the goods, mediation between the seller and the buyer is desirable, but in any case, acceptance of the return remains at the discretion of the Seller and at the buyer's transport costs.

9. Warranties

Seller warrants that all Goods delivered pursuant to an Order shall strictly conform to the terms, conditions, specifications, descriptions, drawings and data specified in the Order or furnished herewith and shall be of good design, material, and workmanship, free from defects, merchantable and fit for their intended purpose and shall meet applicable governmental safety standards.

Seller further warrants that Buyer will have free and clear title to the Goods furnished on the Order and the right to sell such Goods at the time of delivery to Buyer. All Goods shall be new and unused (unless otherwise specified in the Order) at the time of delivery to Buyer. Seller shall also, without any additional charges, assign or otherwise transfer to Buyer all warranties on Goods and services incorporated into the Goods sold on the Order. All warranties herein mentioned shall survive any intermediate or final inspections, delivery, acceptance or payment by Buyer, and all such warranties shall run to Buyer, its successors, assigns, customers, and users of the Goods. No warranties hereunder shall be deemed disclaimed or excluded except in writing signed by an authorized representative of Buyer. If Buyer deems the Goods defective, Seller, at Buyer's election, shall promptly repair or replace the defective Goods, or give Buyer a credit or refund for the price paid for the defective Goods, at any time within Seller's standard European warranty period for the Goods. Alternatively, Buyer may repair defective Goods at Buyer's facility and charge Seller the reasonable cost thereof. Seller shall be responsible for all shipping charges (both to and from Seller's facility) for Goods which are repaired or replaced by Seller.

10. Order Modification

Buyer may at any time, by written notice to Seller, make changes in the specifications, quantities, delivery schedules and shipping instructions in an Order.

Seller will notify Buyer within 10 days if any such change increases or decrease Seller's cost of performing the Order or the time required for its performance in which case an equitable adjustment shall be agreed upon by the Parties. Any stenographic or clerical errors of Buyer in any Order are subject to correction by Buyer.

11. Order Termination

Buyer may at any time terminate for its convenience Seller's performance of an Order, in whole or in part, by written notice to Seller, whereupon Seller shall terminate its performance on such date of notice. If Seller has commenced performance of the Order, Buyer shall negotiate reasonable termination charges with Seller. Seller shall use best efforts to mitigate its losses as the result of any Order termination by Buyer. Buyer may terminate an Order for default with respect to all or any part of undelivered Goods if Seller (a) does not meet the delivery or performance schedule specified in the applicable Order, (b) breached any terms hereof, including warranties, or (c) becomes insolvent or commits or becomes subject to an act of bankruptcy.

12. Set-Off

Buyer may set off any amount Buyer owes to Seller against any amount that Seller owes to Buyer.

13. Infringement

Seller shall indemnify Buyer for any losses, costs, damages and expenses (including court costs and attorneys fees) arising from a claim of infringement of any intellectual property right of a third party relating to any of the Goods delivered. Seller shall defend or settle, at its own expense, any lawsuit or proceeding brought against Buyer by a third party for such infringement. Furthermore, in the event Buyer should be enjoined from the use or sale of the Goods as a result of such lawsuit or proceeding, Seller shall pay all Buyer's costs and expenses including consequential damages and at Buyer's option shall: (a) procure the right for Buyer to continue using and selling the Goods at Seller's expense; (b) replace the Goods with non-infringing Goods at Seller's expense; and/or (c) remove the Goods from Buyer's premises and refund to Buyer the amount paid for the Goods. The provisions of this paragraph shall not apply to any claims, demands, lawsuits or injunctions brought by a third party where the Goods have been manufactured by Seller in accordance with Buyer's specific instructions, specifications, designs or drawings.

14. Indemnification

Seller shall indemnify Buyer from and against any and all losses, costs, damages, and expenses (including court costs and attorneys' fees) arising from any claims against Buyer for (i) injury or death to any individual, or damage to property, alleged to have been caused in whole or in part by defective Goods furnished by Seller or by any act or omission, negligence or otherwise, of Seller or any subcontractor of Seller or of any of their respective employees, workmen, or agents, or (ii) Seller's breach of its warranties and obligations hereunder. Seller will be exempted from any indemnification of the Buyer or any injury, death, or loss caused exclusively by the negligence or willful misconduct of Buyer.

15. Regulatory

Buyer is a medical device manufacturer and is required to comply with different governmental regulatory requirements. If requested by Buyer, Seller agrees to allow Buyer and/or its representatives to conduct an onsite supplier audit of Seller's facilities in connection with Buyer's regulated operations. Alternatively, Buyer may elect to send Seller a mail-in audit which Seller agrees to complete and timely return to Buyer. Seller agrees to cooperate with Buyer in any such audit or any follow up to such audit.

16. Compliance with Laws

In the performance of work hereunder, Seller shall comply with all applicable international, federal, state, and local laws, rules and regulations.

17. Regulation (EU) 2017/745 [MDR]

In case of Medical Devices (MD) the Seller agrees to supply only Goods fully complying with MDR 2017/745 (and its modifications) providing, at Buyer's request, all technical documentation (i.e. EU Declaration of Conformity, Reach / RoHS / Oeko-Tex declarations, instructions, extract of the technical dossier, etc) and all other information useful in order to avoid a misuse by the final Consumer and grant fully compliance to all regulatory duties. Signing this contract the Seller confirms he will supply only Goods manufactured without using raw materials, substances of very high concern (SVHC), semi-finished products and/or other components who may constitute an effective risk/danger (individually or as a whole) according to Regulation (EC) 1907/2006, EN ISO 10993, MDR 745/2017 and their modifications. Every single package must be labeled by the Seller accordingly with MDR 745/2017 using the information, layouts and instructions provided by the Buyer. Non complying Goods will be automatically rejected once the Buyer, during the quality checking of incoming goods, has identified them informing immediately the Seller who'll be in charge for quickly activating the restitution process and replace (at his own expense) the products as soon as possible. In this case, the purchasing order will be consider only partially fulfilled so the Buyer will be entitled to ask for an indemnification for any suffered damage.

18. Assignment

Seller shall not assign any of its rights or interest in this Agreement or any Order or subcontract its performance of this Agreement or any Order, without Buyer's prior written consent. Any prohibited assignment shall be void.

19. Applicable Law

Laws of European Community, without regard to conflict of law principles thereof, will govern this Agreement, and any order issued hereunder.

The Parties expressly consent, and submit themselves, to the exclusive jurisdiction of the Switzerland Confederation for adjudication of any claim.

20. Confidentiality

Seller shall not, except as required in performance of an Order, disclose or use at any time any information which is not publicly known about Buyer's products, processes or services and which is disclosed to or learned by Seller prior to or during such performance and fully comply to the MNDA.

Seller shall, upon demand by Buyer, surrender to Buyer any and all confidential or proprietary documents and other information of Buyer in whatever form which are in Seller's possession or control.

21. Buyer's Information

Buyer retains title to all information and materials (including, but not limited to, all drawings, designs, specifications, technical data and production or product "know-how") furnished to Seller to facilitate performance of an Order, and the same shall be (i) treated as Buyer's confidential information, (ii) used exclusively by Seller to complete the Order, and (iii) returned to Buyer at Buyer's direction or upon completion, termination, or cancellation of the Order, along with all copies or reproductions thereof, unless otherwise agreed to in writing by Buyer.

Seller shall thereafter make no further use of any such drawings, designs, specifications, technical data and production or product "know-how" or of any information derived therefrom in the performance of work for any other customer without Buyer's prior written consent.

22. Seller's Information

Seller retains title to all information and materials (including, but not limited to, all drawings, designs, specifications, technical data and production or product "know-how") furnished to Buyer to facilitate performance of Sales, and the same shall be (i) treated as Seller's confidential information, (ii) used exclusively by Buyer to complete the sales, and (iii) returned to Seller at Seller's direction or upon completion, termination, or cancellation of the Order, along with all copies or reproductions thereof, unless otherwise agreed to in writing by Seller.

Buyer shall thereafter make no further use of any such drawings, designs, specifications, technical data and production or product "know-how" or of any information derived therefrom in the performance of work for any other customer without Seller's prior written consent.

23. No Publicity

Without Buyer's prior written consent, Seller shall not release any news release or publicity regarding this Agreement, an Order, or the fact that Seller is furnishing Goods to Buyer.

24. Non-Waiver

No waiver of any provision of this Agreement or any Order or any right or obligations of either Party shall be effective, except pursuant to a writing signed by the Party waiving compliance, and any such waiver shall be effective only for the specific instance and purpose stated in such writing

25. Final provisions

For all that is not regulated in this document, the European regulations in force at the time of conclusion of the contract are valid.



ORTHOSERVICE RO+TEN

Take care feel better



www.orthoservice.com



UNI CEI EN ISO 14001:2021
UNI EN ISO 9001:2015



9 DCCT00048 INT01.0323

Headquarter: ORTHOSERVICE AG
Via Milano 7 · CH-6830 Chiasso (TI) · Switzerland
Tel. 0041 (0) 91 822 00 88 · Fax 0041 (0) 91 822 00 89
info@orthoservice.com · www.orthoservice.com
Niederlassung Deutschland: Orthoservice Deutschland GmbH
Flugstraße 8 · D-76532 Baden-Baden · Deutschland
Tel. 0049 (0) 7221 991 39 11 · Fax 0049 (0) 7221 991 39 13
info@orthoservice.de · www.orthoservice.de
Siedziba w Polsce: ORTHOSERVICE POLSKA Sp. z o.o.
ul. Warszawska 416a · 42-209 Częstochowa Polska
Tel.: +48 (0) 34 340 13 10 · (NIP): 9492246785
www.orthoservice.pl · info@orthoservice.pl
Sede italiana: RO+TEN S.r.l.
Sede legale: Via Marco De Marchi, 7 · I-20121 Milano (MI) · Italia
Sede operativa e amministrativa:
Via Comasina, 111 · I-20843 Verano Brianza (MB) · Italia
Tel. 0039 039 601 40 94 · Fax 0039 039 601 42 34
info@roplusten.com · www.roplusten.com

Società soggetta a Direzione e Coordinamento (art. 2497bis CC):
Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland

Copyright Orthoservice Ro+Ten March 2023

The descriptions and images in this document are for illustrative and commercial purposes only. Orthoservice reserves the right to make modifications without prior notice based on their needs.

ORTHOSERVICE AG

Via Milano, 7
CH-6830 Chiasso (TI)
Switzerland
Tel. +41 (0) 91 822 00 88
fax +41 (0) 91 822 00 89
info@orthoservice.com
www.orthoservice.com

